

ZerofloX™

THE NEW STANDARD FLOCK-FREE
MICRO APPLICATOR

ZERO IS
HERO



 medmix

DENTAL 



If you are looking for **ZERO compromises – look no further**

The ZerofloX™ micro applicator sets a new standard compared to conventional micro applicators. Its flock-free, pliable brush head and ergonomic handle with two bending points deliver precision application for esthetically pleasing clinical results.

Say hello to ZerofloX. The new standard flock-free micro applicator.



Soft, pliable head
Fiber-free
elastomer bristles




**Two bending
points**



Grip handle



**Triangular-shaped
handle**




**Customized
colors**

1 ZERO FLOCKS

ZerofloX™ is FLOCK-FREE. Its fiber-free elastomer bristles eliminate loss of adhesively-fixed fiber flocks, which could remain in the final restoration as debris*.

The results do not simply look good – they look genuinely perfect.

* Journal of Dentistry, Vol. 127, Dec. 2022, 104320
<https://doi.org/10.1016/j.jdent.2022.104320>




Clinical photo showing lost fiber flocks using a conventional flocked micro applicator.
Courtesy of Prof. Fabio M. Salomão DDS, MsC, PhD, Brasil.

3 ZERO LOSS

ZerofloX is EFFICIENT. Its non-absorbing, drip-free elastomer bristles uptake and release a consistent amount of material without depositing excess.

With its optimized efficiency, ZerofloX minimizes material waste and maximizes cost savings with every application.



Data available upon request.

2 ZERO TROUBLE

ZerofloX is PRECISE. Its soft, pliable head allows easy intraoral access in difficult-to-reach, tight-fitting or irregular areas. In addition, the triangular-shaped handle with two bending points promises ergonomic use.

Treating patients is easy and precise.



Courtesy of Prof. Fabio M. Salomão DDS, MsC, PhD, Brasil.

4 ZERO CONFUSION

ZerofloX is UNIVERSAL. Its flexible but dimensionally stable, single-diameter head size can be used in combination with a variety of clinical applications and dental materials.

A true one-size-fits-all solution that makes application as straightforward as possible!



ZerofloX micro applicator positioned inside the root canal of an extracted tooth.
Courtesy of Prof. Fabio M. Salomão DDS, MsC, PhD, Brasil.

ZERO IS HERO

Save yourself from compromising on quality, comfort and efficiency. ZerofloX™ is your chance to leave second best behind.

Please contact us to get your free ZerofloX™ sample today.

What dentists are saying



"With its fiber-free monobloc elastomer bristles, the ZerofloX micro applicator has the potential to become the future gold standard in adhesive esthetic dentistry."

Prof. Dr. med. dent. Ivo Krejci,
Switzerland



"ZerofloX sets itself apart from all existing micro applicators and will eventually represent a new standard for the category. During clinical use, my patients experience greater comfort when applying intraoral products. ZerofloX provides better surface wettability and superior durability, allowing repeat use in more than one tooth or restoration, eliminating the need to use more than one micro applicator per patient. This drastically reduces the number of disposable applicators in the clinical routine."

Prof. Paulo Vinicius Soares, DDS, MS, PhD,
Brazil



"The ZerofloX micro applicator offers significant advantages compared to conventional flocked-type applicators due to its ergonomic brush head design, which allows for easy access and application onto just about every type of tooth surface/cavity without the loss of hairlike fiber flocking."

Prof. Fabio M. Salomão DDS, MSc, PhD
Universidade Estadual de Londrina
São Paulo, Brasil

medmix Switzerland AG
Ruetistrasse 7
9469 Haag, Switzerland
P +41 81 414 70 00
www.medmix.swiss



Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC, ZerofloX and Candy Colors are trademarks of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland AG 2023